

The Encode® Emergence Impression System


Healing Abutments

Surgical Placement


The [Encode Emergence Healing Abutment](#) is a 3-in-1 two-piece healing abutment/impression coping/scan-body designed to facilitate gingival tissue healing and consists of an abutment and retaining screw that are assembled together.

Encode Emergence Healing Abutments have the Encode coding scheme on the occlusal surface and a hex connection at the base of the healing abutment that engages the hex with the implant for orientation and anti-rotation.


Grooves on the occlusal surface of the screw identify the implant system.



TSV
3 Grooves



Certain
2 Grooves



Eztetic
1 Groove

1



Following implant placement, select an Encode Emergence Healing Abutment with the appropriate restorative platform diameter, emergence profile, and collar height.

2



Use wide Encode Emergence Healing Abutment(s) to ensure the final abutment will seat easily without excessive blanching of the soft tissue. Bone profiling may be required to fully seat the healing abutment, particularly in cases with sub-crestal implant placement deeper than 1 mm.



Note: Please ensure that the instructions for use are made available to the restorative clinician and use the patient chart labels located in the package of the Encode Emergence Healing Abutment.



Use an Encode Emergence Healing Abutment that stands taller than the surrounding soft tissue so all sides of the occlusal Healing Abutment are clear to facilitate proper traditional or digital impression and subsequent scanning.

Implant System	Driver Part Number
TSV/Eztetic	HXGR1.25 and HXLGR1.25
Certain	PHD02N and PHD03N

Pick up the Encode Emergence Healing Abutment(s) using a hexed driver (see list in table above). When placing the Encode, please ensure that the hex is oriented and engaged prior to driving the screw all the way in for proper seating. Place the healing abutment on the implant and hand-tighten.



Following placement, confirm that the Encode Emergence Healing Abutment is fully seated on the implant with a radiograph.



Once fully seated, torque the healing abutment screw to 20 Ncm using a restorative torque wrench and torque wrench hex driver (see list in table below). Where applicable, suture the soft tissue around the Encode Emergence Healing Abutment.

Implant System	Torque Wrench	Driver Tip
TSV/Eztetic	TWR	TW1.25 or TW1.25L
Certain	L-TIRW	IHX1.2D, RASH3N, IXL1.2D or RASH8N

For more information visit ZimVie.com

ZimVie US Corp LLC
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272
dentalCS@ZimVie.com



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. and its affiliates; and all products are manufactured by one or more of the dental subsidiaries of ZimVie US Corp LLC (Biomet 3i, LLC, etc.) and distributed by ZimVie US Corp LLC and its authorized marketing partners. EU Medical device classification: Dental healing abutment (Class IIb). The Intended purpose of a healing abutment is to shape soft tissue during implant healing. The products here shown are not eligible for reimbursement in France. For other countries, please check the local laws and regulation. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZVINST0005-GLBL REV D 11/25 ©2025 ZimVie. All rights reserved.

