




Stability and confidence
in your digital workflow



Implant Analogs for TSX[®], TSV[®], and Eztetic[®] Implant connections

Innovative Design

Each Implant Analog utilizes self-locking and anti-rotational features for greater stability. When used in 3D-printed models, a single fixing 1.27 mm hex screw is included, ensuring the analog is securely attached allowing accurate positioning.




Precise Matching Connection

TSV Implant Analogs are color-coded to match the platform and are manufactured in titanium alloy. Analogs can be placed in traditional stone and 3D models using the same tooling that you are familiar with for TSV.



Implant Analogs for TSX and Eztetic Implant Systems




Color-coded by implant platform

Implant Platform	Item No.
 NP (2.9 mmD)	NPILA-29
 3.5 mmD	ILATSV35
 4.5 mmD	ILATSV45



Implant Analogs for TSV Implant Systems

Color-coded by implant platform

Implant Platform	Item No.
 3.5 mmD	ILATSV35
 4.5 mmD	ILATSV45
 5.7 mmD	ILATSV57



For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272



Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental Inc., etc.) and marketed and distributed by ZimVie Dental (formerly Zimmer Biomet Dental) and its authorized marketing partners. TSV Implant Analogs are manufactured by Terrats Medical S.L. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians use only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1284 REV B 04/24 ©2024 ZimVie. All rights reserved.

