



Protect. Heal.  
Repair.



# Collagen Matrices

## Plug, Tape, and Patch



Absorbable wound dressings

### 1 Protects Wound Bed

Adheres and provides coverage to oral wounds and sores

### 2 Aids Wound Healing

Porous, absorbable matrix supports delicate new tissue

### 3 Highly Porous, Highly Absorbent

ZimVie's Collagen Plug holds up to 30x own weight in fluid  
Resorbs in fewer than 30 days  
Greater than 90% open pores

### 4 Versatile for Everyday Use

Three convenient shapes for common oral wounds and procedures



# Absorbable Collagen Matrices

In three convenient shapes

Collagen Plug



10 mm x 20 mm  
Cat. No. 0102Z

Collagen Tape



2.5 cm x 7.5 cm, 1 mm thick  
Cat. No. 0100Z

Collagen Patch



2 cm x 4 cm, 3 mm thick  
Cat. No. 0101Z

## Versatile wound protection

- Designed to resorb within a short timeframe
- Indicated for management of oral wounds and sores
  - Denture sores
  - Oral ulcers (non-infected or viral)
  - Periodontal surgical wounds
  - Suture sites
  - Burns
  - Extraction sites
  - Surgical wounds
  - Traumatic wounds



Atraumatic tooth extraction



Modified Collagen Plug delivered to protect grafted site



Sutured site



Ovate pontic contours soft tissue during healing

For more information, visit [ZimVie.com](http://ZimVie.com)

ZimVie  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454

 **ZimVie**  
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Zimmer Collagen Plug, Tape and Patch are manufactured by Collagen Matrix, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1919 REV A 11/23 ©2023 ZimVie. All rights reserved.

